DS-1031

**NISSEI**<sup>®</sup>

JAPAN

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
This manual is intended to assist you in the safe and efficient operation of BLOOD PRESSURE MONITOR DS-1031. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

## PART NAMES AND PRODUCT COMPONENTS

1. Cuff
2. Power Adapter
3. Batteries
4. Monitor

- A. Air plug
- B. Battery Compartment
- C. AC Connector
- D. "Set" Key
- E. Memory Key 1
- F. Memory Key 2
- G. "Start/Stop" Key
- H. Display
- I. Air connector
- K. Air hose

## GENERAL INFORMATION

### Indications for use

This product is intended for noninvasive measurement of systolic and diastolic blood pressure, determination of pulse rate and calculation of pulse pressure in adults in a home healthcare environment. The product is not designed for neonatal use. Please consult with your doctor or physician to use this product to take blood pressure of child or person in pregnancy or under pre-eclamptic condition.

The cuff is applicable to arm with circumference between 22 and 42 cm. Pressure is measured in the range from 40 to 250 mm Hg, and pulse rate in the range from 40 to 160 beats per minute.

## Method Of Measurement

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The cuff is connected to the main unit and wrapped around the arm. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value.

## Used new NISSEI Technologies

**Fuzzy Inflation** – This product employs FUZZY inflation to ease stress on wrist during cuff inflation. By detecting pulse signals while inflating, this product is able to determine the sufficient pressure for each measurement.

**Arrhythmia Indicator** – displays when the pulse rhythm detected during measurement was irregular. Pulse rhythm can be disturbed from talking, moving or arrhythmias.

**Touch Control** – touch-sensor operation.

**Body Motion Indicator** – blood pressure value taken while moving cannot be said to be the correct value because body movement can affect blood pressure.

**Cuff Symbol** – displays whether the cuff is applied properly or not.

**Reliability Symbol** – displays when the measurement is supposedly made in good condition.

**Pulse Pressure Display** – calculates and displays pulse pressure. Pulse pressure is value subtract diastolic value from systolic value and therefore pulse pressure tends to increase with age.

**ATTENTION!** Do not use cuffs other than the original cuff included with this product.

## COMPLETE SET

1. Main body – 1 pcs.
2. Cuff – 1 pcs.
3. Battery – 4 pcs.
4. AC Adapter – 1 pcs.
5. Carrying bag – 1 pcs.
6. Instruction manual – 1 pcs.
7. Warranty card – 1 pcs.
8. Packing – 1 pcs.

## TIPS ON TAKING YOUR BLOOD PRESSURE

- 1 Do not use this product without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.
- 2 Use of this product in areas near mobile phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.
- 3 For correct measurement it is necessary to know that ARTERIAL BLOOD PRESSURE SUBJECTED TO SHARP FLUCTUATIONS IN SHORT GAP OF TIME. The level of the arterial blood pressure depends on many factors. Usually level of arterial blood pressure is lower in summer and higher in winter. Arterial blood

pressure changes together with atmospheric pressure and depends on many factors, e.g. physical loads, emotional excitability, stress, meals, dietary pattern, etc. Different medicines, alcohol, smoking render large influence to level of arterial blood pressure. When blood pressure is measured in a hospital, it could be higher than when measured at home. This is because you are tense at the hospital and relaxed at home. Because arterial blood pressure increases in low temperature, determine arterial blood pressure in room temperature (approximately 20° C). If this device was stored in low temperature, it is necessary to bear it before using at least 1 hour in room temperature; otherwise, result of the measurement can be wrong. Blood pressure changes easily and the difference in persons health during day can be before 30-50 mmHg of systolic (higher) pressure and before 10 mmHg of diastolic (lower) pressure. The dependency of the arterial blood pressure on different factors is individual in every person. Therefore, it is recommended to write the special diary book of data of arterial blood pressure. **ONLY QUALIFIED PHYSICIAN CAN INTERPRET YOUR ARTERIAL BLOOD PRESSURE DATA OF YOUR DIARY BOOK.**

4 Make arterial blood pressure measurement in time determined by your personal doctor in case of cardio-vascular diseases and other diseases when monitoring of arterial blood pressure is necessary. **REMEMBER! DIAGNOSTICS AND TREATMENT OF HYPERTENSIONS CAN BE MADE ONLY BY QUALIFIED PHYSICIAN ON BASE OF ARTERIAL PRESSURE DATA GOT BY PHYSICIAN ONESELF. ACCORDING ADMINISTRATION OR CHANGE OF MEDICINE DOSAGES IT IS NECESSARY TO FOLLOW TO PRESCRIPTION MADE ONLY BY YOUR DOCTOR.**


Fig.1

5 The right measurement of arterial pressure can be difficult in deep arteriosclerosis, in weak pulse wave, in patients with expressed changes of heart contraction rhythm. Let a qualified physician interpret your blood pressure readings.

6 **IT IS NECESSARY TO KEEP THE SILENCE DURING MEASUREMENT BY THIS DEVICE TO GET CORRECT DATA OF YOUR ARTERIAL BLOOD PRESSURE.** Measurement of the arterial blood pressure should be conducted in a quiet environment in room temperature. Avoid eating for 1 hour, drinking alcohol for 1,5-2 hours, smoking and other activities that affect your blood pressure right before a measurement.

7 For accuracy of the arterial blood pressure measurement the measuring range of the cuff should be corresponded to the size of your upper arm. **THE CUFF SHOULD NOT BE SMALLER OR ON THE CONTRARY GREATER THAN NECESSARY.**

8 Make repeated arterial blood pressure measurements with interval 3 minutes to restore the blood circulation. Suffering from pronounced atherosclerosis persons are required more time interval of the measurements (10-15 minutes) because elasticity of vessels decreased significantly in this disease. That is same for patients suffering from sugar diabetes for a long period of time. For more exact determination of the arterial pressure, conduct the series of 3 consequent measurements then to calculate the average of measurement results.

## BATTERY INSTALLATION


Fig.2


Fig.3

1. Open the battery compartment cover (fig.2).
2. Install four "AA" type batteries into the compartment. Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment (fig.3). Batteries can be easily installed or removed as their (-) ends are pushed against the spring.

You may use nickel hydride rechargeable batteries. These batteries can not be charged with this product. Use the designated battery charger.

3. Close the battery compartment cover.

Do not force the battery cover into position.

**Battery symbol** will flash when battery power is low.

Measurement can not be made while the symbol is displayed. Replace the batteries. All four batteries need to be replaced with new ones; do not mix new and old batteries or different batteries.

The symbol might appear only during measurement although you can review readings in memory. This is because more battery power is required for conducting a measurement, e.g. inflating the cuff, than displaying the readings in memory. Please have the batteries replaced. The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.

The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly.

## Operation With AC Adaptor

AC Adapter socket is located at the side of the device. Connect the AC Adapter jack to the device firstly and after supply AC Adapter plug into the power outlet, press «START / STOP». To switch off the device after finishing the measurement press «START / STOP», unplug the AC Adapter from the power outlet and disconnect AC Adapter from the device.

### WARNING!

If you will switch off AC Adapter without batteries in the device all information stored in memory (results of measurements, date and time) will be erased. If you want to save this information, do not remove batteries from the device when using the AC Adapter.

## ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

Date and time can be set after installing the batteries.

Set the time guarantees the preservation of the measurement results with the correct date and time. There is a possibility of measurement of pressure and pulse rate without setting date and time.

Touch "SET" KEY and keep your finger touched to the key until "2011" starts flashing.

Clock is set in the order of year, month, day, hour and minute.

The flashing number increases with  and decreases with . The number will be fast-forwarded if you keep your finger on the key. Touching "SET" will fix the number and the next article will flash. Touching "START/STOP" will terminate the setting.

**IMPORTANT!** If the date and time have been set, current time will be shown on the display when the device is turned off.

## CORRECT MEASURING POSTURE

Sit on a chair. Lift your left hand slightly palm up and place your elbow on the table. Place the cuff at heart level by placing the forearm case or a folded towel (fig.4).

Pressure measurement while lying. Place the cuff at heart level, using a case or a folded towel (fig.6).


Fig.4

Fig.5


Fig.6

The results of measurements can vary slightly depending on the posture during the measurement. If the cuff is below (above) in relation to the heart, the measurements will be overstated (understated).

## WRAPPING THE CUFF


Fig.7


Fig.8


Fig.9


Fig.10


Fig.11

1 Put the cuff on your left arm with AIR HOSE positioned toward your hand. (fig.7). You may use your right upper arm for measurement (the results of measurements can slightly).

2 Attach the cuff with hook and loop fastener with the lower edge of the cuff approximately 2 to 3 cm above the inner elbow. Adjust the position of the cuff so that AIR HOSE is over the inner part of your arm over the brachial artery (fig.8).

3 Press the surface of the cuff to make sure that the hook & loop fastens securely (fig.9). If the cuff is wrapped tightly or loosely, inaccurate blood pressure readings may result.

4 If the arm has taper form, it is recommended to put the cuff on "spiral way", as shown in the figure (fig. 10).

5 If you are wearing a shirt that might restrict circulation in your upper arm or you roll your sleeve up over the upper arm, the blood flow will be restricted, preventing accurate measurement (fig.11).

## PROCEDURE FOR THE MEASUREMENT

**IMPORTANT!** The device has sensor keys and clicking is a slight touch of a finger.

Moisture, dirt and other objects between the finger and panel of the device can affect the ability to respond to the touching keys.

1. Insert AIR PLUG into AIR CONNECTOR before starting a blood pressure measurement.

*Take deep breaths and relax. Do not move, chat or strain your arm or hand during measurement.*

2. Touch "START/STOP". Air is exhausted from the cuff (fig.12).

3. Automatic inflation starts (fig.13).

4. Inflation stops when the pressure, that is, the displayed value, reaches the designated value.


Fig.12


Fig.13



Fig.14



### Body Motion Symbol


Blood pressure value taken while moving cannot be said to be the correct value because body movement can affect blood pressure. This product analyzes pulse wave and displays «» when body motion is detected. «» indicates the results might be affected by body movement.

### Cuff Symbol

This product displays whether the cuff is applied properly or not by analyzing the pressure wave. «» indicates that the cuff is applied properly but «» indicates that the cuff may be wrapped too tightly or too loosely.

Touch "START/STOP" again to cancel measurement. The monitor will exhaust air from the cuff and turn off.

5. Pulse Symbol «» is displayed as pulse is detected (fig.14).

6. Air is released from the cuff and measured values are displayed as measurement completes (fig.15).

7. Touching  or will  switch the bank.

8. Touch "START/STOP".

The monitor will be turned off. Even if you do not turn off the monitor, it turns off automatically after 3 minutes.



Fig.15


### Reliability Symbol

Some of the important aspects in measuring blood pressure are how properly the cuff is applied and how calm you are during measurement. These aspects can affect measurement values and such values cannot be said to be correct values. This product shows «» when the measurement is supposedly made in good condition without influences from such aspects.

Do not execute repeated measurements for congestion of blood could result in false measurement. Let your wrist rest for at least 5 minutes.

## Irregular Pulse Rhythm Symbol

Pulse rhythm can be disturbed from talking, moving or arrhythmias. This product displays «⊕» when the pulse rhythm detected during measurement was irregular (fig.16).

## WHO Classification Symbol

Measured blood pressure values are classified against WHO guideline. Scale, being lower right on display, allows us to estimate the obtained figures according to the classification (fig.17).

Fig.17

WHO scale is the scale classification derived values of blood pressure

Display	WHO Classification	SYS	DIA
	<b>Hypertension (severe)</b>	<b>≥180</b>	<b>≥110</b>
	<b>Hypertension (moderate)</b>	<b>160-179</b>	<b>100-109</b>
	<b>Hypertension (mild)</b>	<b>140-159</b>	<b>90-99</b>
	<b>High Normal</b>	<b>130-139</b>	<b>85-89</b>
	<b>Normal</b>	<b>120-129</b>	<b>80-84</b>
	<b>Optimal</b>	<b>&lt;120</b>	<b>&lt;80</b>

## Pulse Pressure Display

This product calculates and displays pulse pressure (fig.18).

While systolic blood pressure keeps increasing with age, diastolic blood pressure tends to start decreasing at around age 50. Pulse pressure is value subtract diastolic value from systolic value and therefore pulse pressure tends to increase with age. Pulse pressure has been considered to be related to stiffness of artery and has been studied as one of circulatory risk factors. It is sometimes said that 45 mmHg of pulse pressure is normal value. Pulse pressure only does not represent arterial sclerosis. However, observation of pulse pressure in long term will be important.


Fig.18

## MEMORY FUNCTION

The measured values are automatically saved for later review in either one of two memory banks. These two memory banks can be used to save readings of two persons separately or to save morning readings and evening readings separately.

Each bank can save up to 60 readings. When the number of stored readings reaches 60, the oldest reading will be deleted to record new reading.

Saved readings are recalled with measurement date and time when clock function is activated. \*The clock needs to be activated before measurement for date and time to be saved with the measured value.

[Err] results are not saved.

1 Turn off the monitor once with "START/STOP" when the results are displayed after a measurement. Touch [1] to review readings saved in memory 1 and [2] to review those in memory 2.

The average of the saved readings are displayed (fig. 19).

Average will not be displayed unless there are two or more readings saved.

2 Touch the same memory key again. The latest reading saved in the selected bank is displayed.

3 The display at the top shows memory number, measurement date and time alternately.

4 The latest reading is displayed as memory "No. 1". The bigger memory number indicates the older number.

Even if you do not turn off the monitor, it displays the result for approximately 30 seconds and will turn off automatically.

Touching [2] at display of readings in the bank M1 will switch the display to reading in bank M2 and touching [1] side again will return to the display of result in the bank M1.

5 Touch "START/STOP".

### Deleting saved readings

The readings can be deleted from the memory. The memory can be initialized only when two or more results are saved in the bank.

1 Touch [1] to delete reading in MEMORY 1 or [2] to delete reading in MEMORY 2. Display the average reading of to delete entire readings in the memory.

2 The display will start flashing. Keep your finger on the key until the reading is replaced with «---».




Fig.19

## TROUBLESHOOTING

SYMPTOM	PROBABLE CAUSE	Remedies
Blood pressure is extremely high or low.	The cuff is not adjusted to your heart. The cuff is not wrapped appropriately. Movement or speaking was detected.	Adjust the cuff to the height of your heart. Reconfirm how to put the cuff on your arm. Do not move or speak during measurement.
Measured values vary all the time.	Mental and physical condition or measurement condition is influencing blood pressure.	Make measurement under the same condition.
Measured value is different from that taken at a hospital.	Mental condition such as nervousness at a hospital or relaxed feeling at home is influencing blood pressure.	Take records of blood pressure taken at home and consult with your doctor.
Inflation is repeated.	The cuff is inflated again when the initial pressure was insufficient for blood pressure to be taken or movement was detected. Inflation is repeated until blood pressure is taken.	Repeated inflation is not a malfunction of the monitor. Do not move or speak during measurement.
	Over pressure: blood pressure could not be taken due to moving or speaking although the cuff was inflated to the maximum pressure.	Do not move or speak during measurement.
	Blood pressure could not be taken due to moving or speaking.	Do not move or speak during measurement.
	The cuff is not wrapped appropriately.	Reconfirm how to put the cuff on your arm.
	Batteries are weak.	Replace all batteries with new ones.
Clock is not displayed.	Batteries are deleted. Batteries are inserted incorrectly. Battery terminals are not clean. Touching the keys with tip of finger or nail is touching the key.	Replace all batteries with new ones. Reinsert the batteries correctly. Clean the terminals with dry cloth. Touch the key with the plane of your finger.
Measurement date and time are displayed with «--/--».	The clock is not activated. Measurement was taken before the clock was activated.	Adjust date and time and activate the clock. Measurement date and time can not be saved without clock being activated.

Nothing is displayed.	Batteries are deleted. Batteries are inserted incorrectly. Battery terminals are not clean. Touching the keys with tip of finger or nail is touching the key.	Replace all batteries with new ones. Reinsert the batteries correctly. Clean the terminals with dry cloth. Touch the key with the plane of your finger.
	START/STOP was being touched at installation of batteries.	Turn off the monitor once with "START/STOP" and make measurement.

If you cannot get correct measurement with the methods above, contact your dealer. Do not disassemble or tamper with the internal mechanism.

## WARRANTY

- 1 Warranty period for the electronic unit is 5 years from the date of sale. The warranty period for the cuff and power adapter is 12 months from the date of sale.
- 2 The warranty obligations are prescribed by warranty certificate for buyer. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual.
- 3 The addresses of organizations for guarantee maintenance are present in the warranty certificate.

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## SPECIFICATIONS

Operating Principle	Oscillometric
Indicator	15 digits liquid crystal display
Pressure Indicating Range: cuff pressure, mmHg	0-300
Measuring Range: cuff pressure, mmHg. pulse rate, bpm.	40-250 40-160
Accuracy: cuff pressure, mmHg. pulse rate, %	±3 ±5
Inflation	Automatic inflation (Fuzzy Inflation)
Deflation	Automatic (electric control valve)
Power supply	6V, 4 x AA (dry alkaline elements) or AC Adapter
Power Consumption, W	4
Memory	2 x (60 + average)
<b>AC Adapter ADP-W5</b>	
Output voltage, V	6
Maximum electric load, A	0,5
Input voltage range, V/Hz	100-240/50
Operating Environment temperature, °C relative humidity, % Rh	from 10 to 40 85 or below
Environment temperature, °C relative humidity, % Rh	from -20 to 50 85 or below
Cuff Model	Cuff DS-1031
Cuff Size	Universal Slim Fit Cuff (Arm Circumference 22-42 cm)
Overall dimensions: Size (without cuff), mm Weight (without package, case, battery and power adapter), g	115 x 115 x 67 406
The service life: Unit (without cuffs), years Cuff, years	7 3
Year of manufacture:	Year of manufacture is indicated in the serial number after the letters «SN» on the case of the device
Water resistance	IPX0 (IEC 60601-1)
Key to symbols	<p> Type BF equipment</p> <p> Refer to instruction manual/booklet</p> <p> The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly.</p>

## CARE, STORAGE, REPAIR AND RECYCLING

- 1 This appliance must be protected from excessive moisture, direct sunlight, shock and vibration. **UNIT IS NOT WATERPROOF!**
- 2 Do not store or use the device in close proximity to heaters and open flames.
3. When the product has been stored at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using.
- 3 If the product is not used for a long time, remove the battery from the unit. The flow of electrolyte from batteries can cause product's damage. **KEEP BATTERIES OUT OF REACH OF CHILDREN!**
- 4 Do not pollute the unit and keep it dust. For cleaning, you can use a soft, dry cloth.
- 5 Do not let the unit and its parts with water, solvents, alcohol, gasoline.
- 6 Keep the cuff, bulb and the tubes from sharp objects.
- 7 Do not expose the unit to strong shocks and do not drop it.
- 8 If you need to repair only in specialized organizations.
- 9 On the expiry of the period of service should periodically consult specialists (specialized repair organizations) to check the technical condition of the unit.
- 10 When disposing guided by the rules in your area. Special conditions for the disposal of this product is not established by the manufacturer.
- 11 The cuff is resistant to multiple sanitization. Processing is permitted inside the tissue covering the cuff (in contact with the patient's arm) with a cotton swab dipped in 3% hydrogen peroxide solution. With prolonged use, may be a partial bleaching of textiles to the cuff. Do not wash the cuff, as well as treatment with a hot iron.
- 13 Do not leave the device unattended when it plugged into the power outlet.

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© Copyright 2011-2012.
# NIHON SEIMITSU SOKKI CO., LTD.
JPI Inc. EU Office Neubertstrasse 32,